



Acetabular Defect Reconstruction In Hip Revision

PETER BREHM Revision Cup Algorithm


Classification acc. to Paprosky




Type I: No significant osteolyses, no implant migration, preserved acetabular rim, columns intact




Type IIa: Cranomedial migration of the implant, cavitory defect, columns intact




Type IIb: Cranomedial migration of the implant, segmental defect, columns intact



Type IIc: Medial migration of the implant, medial defect, moderate osteolysis, impaired anterior column



Type IIIa: Severe migration, anterior column intact, posterior column impaired



Type IIIb: Severe migration, anterior and posterior column impaired, pelvic discontinuity