

BPK-S INTEGRATION REVISION

SIZING INTEROPERABILITY

Sizing interoperability Femur SC, Tibia, Patella and Insert

FEMUR SIZE = PATELLA SIZE = INSERT SIZE	TIBIA SIZE
Sz. 3	Sz. 3, 4, 5
Sz. 4	Sz. 3, 4, 5, 6
Sz. 5	Sz. 4, 5, 6, 7 ¹
Sz. 6	Sz. 6, 7 ¹ , 8 ¹

¹ Sizes 7 and 8 of Tibia are only available in versions UC/SC and SC.

Sizing interoperability Femur RH/TH, Tibia, Patella and Insert

FEMUR SIZE = PATELLA SIZE = INSERT SIZE	TIBIA SIZE
Sz. 3	Sz. 3, 4, 5
Sz. 4	Sz. 3, 4, 5, 6
Sz. 5	Sz. 4, 5, 6
Sz. 6	Sz. 6

The size of the patella and the insert comply with the size of the femur. For example a femur size 4 has to be used with a patella and insert component of the size 4.

Patella components available

FEMUR SIZE = PATELLA SIZE	PATELLA DIAMETER
Sz. 3	Ø 24 mm Ø 28 mm Ø 32 mm
Sz. 4	Ø 24 mm Ø 28 mm Ø 32 mm Ø 36 mm
Sz. 5	Ø 28 mm Ø 32 mm Ø 36 mm
Sz. 6	Ø 28 mm Ø 32 mm Ø 36 mm

Patella height is 8 mm for all diameters.