

BPK-S INTEGRATION

SIZING INTEROPERABILITY

Sizing interoperability Femur, Tibia, Patella and Inlay

FEMUR SIZE = PATELLA SIZE = INLAY SIZE	TIBIA SIZE
Sz. 1	Sz. 1, 2, 3
Sz. 2	Sz. 1, 2, 3, 4
Sz. 3	Sz. 2, 3, 4, 5
Sz. 4	Sz. 3, 4, 5, 6
Sz. 5	Sz. 4, 5, 6, 7
Sz. 6	Sz. 6, 7, 8
Sz. 7	Sz. 6, 7, 8

The size of the patella and the inlay comply with the size of the femur. For example a femur size 4 has to be used with a patella and inlay component of the size 4.

Patella components available

FEMUR SIZE = PATELLA SIZE	PATELLA DIAMETER
Sz. 1	Ø 24 mm Ø 28 mm
Sz. 2	Ø 24 mm Ø 28 mm
Sz. 3	Ø 24 mm Ø 28 mm Ø 32 mm
Sz. 4	Ø 24 mm Ø 28 mm Ø 32 mm Ø 36 mm
Sz. 5	Ø 28 mm Ø 32 mm Ø 36 mm
Sz. 6	Ø 28 mm Ø 32 mm Ø 36 mm
Sz. 7	Ø 28 mm Ø 32 mm Ø 36 mm

Patella height is 8mm for all diameters.